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Marca|Trademark:Nahita

The new Nahita paraffin dispensing unit is made up of an intuitive control panel, an upper paraffin reservoir and dispenser, mould and cassette warmers, as well as a large working area. With an ergonomic and compact design, this equipment is ideal and essential in histology and pathological anatomy laboratories for embedding tissues in paraffin blocks for subsequently being cut with the microtome.

#### Control panel

Clear and easy to understand, it allows an independent temperature adjustment of the different regions of the equipment (working area, paraffin reservoir, mould and cassette warmers and dispenser). Temperature of each one of these parts can be adjusted between 55 °C and 70 °C in steps of 5 °C. Furthermore, the control panel presents an easy-to-read LCD display and allows programming the operation start at the desired hour and week day so as to have always the equipment ready to use.

#### Working area

The working platform is situated at an ideal height to avoid fatigue after long periods of work and presents a large easy to clean granite surface, a region with independent temperature adjustment and a forceps holder.

The heated region is provided with a drainage system to accumulate the left-over paraffin into a collecting tray placed just under the working platform and easily removable for emptying. The heated forceps holder, with capacity for 6 forceps, is easily accessible from both sides.

For a higher user's comfort, the working area can be illuminated by a non blinding light to facilitate sample orientation inside the paraffin block.

#### Tanks

Paraffin reservoir is located at the upper part of the equipment and presents a capacity of 2 L, minimizing the need for constant refilling. The equipment also presents two tanks to keep moulds and tissue cassettes warm at any time.

Each tank is provided with a metallic lid to keep liquid paraffin free of dust and particles.

#### Paraffin dispenser

According to user necessities or preferences, the dispenser can be activated either manually or via the foot switch supplied with the equipment. The manual dispenser clip is smooth-running to avoid fatigue after long periods of work.