# AGITADOR DE VARILLAS DIGITAL DIGITAL OVERHEAD STIRRER AGITATEUR DIGITAL À TIGE

Ref. | Code | Réf. LLB008 - Mod. 688

The logo for Nahita Blue, featuring the word "nahita" in a bold, blue, sans-serif font above the word "blue" in a smaller, lighter blue, sans-serif font.

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*This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.*

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**1. INTRODUCTION**

The Nahita-Blue Model 688 overhead stirrer features an advanced international design and state-of-the-art manufacturing technology. It has a compact and robust structure, with a modern and functional design. Its high versatility and ease of use make it an ideal device for various types of scientific applications.

This stirrer is specifically designed for mixing high-viscosity liquids, such as liquid mixtures and solid-liquid mixtures.

**2. TECHNICAL SPECIFICATIONS**

<b>Code</b>	LLB008
<b>Model</b>	688
<b>Control</b>	Digital display with buttons
<b>Speed range</b>	100-1200 rpm
<b>Speed increment</b>	10 rpm
<b>Timing</b>	0-24 h
<b>Maximum capacity</b>	5 l
<b>Maximum viscosity</b>	5000 mPas
<b>Motor</b>	Brushless motor
<b>Maximum torque</b>	10 N.cm
<b>Chuck diameter range</b>	1-10 mm
<b>Power</b>	36 W
<b>Voltage</b>	100-220V 50/60 Hz
<b>Dimensions</b>	134x70x173 mm

### 3. OPERATION

- **Application:** It is mainly used for stirring high-viscosity liquids.
- **Usage:** It is widely used as an experimental tool in higher education institutions, environmental protection, health, epidemic prevention, the chemical industry, tap water facilities, and other applications.
- **Using steps:**
  - Connect the device and switch it on
  - Press 'MENU' to set the speed and time parameters
  - Adjust the desired parameters using the up and down buttons
  - Press 'START' to begin
  - Press 'STOP' to stop operation
- **Instrument use environment:**
  - Indoor use only
  - Altitude:  $\leq 2000$  m
  - Working temperature range:  $+5$  °C a  $+40$  °C
  - Humidity range:  $\leq 80\%$
  - No vibrations or airflow that could affect performance.
  - No conductive dust, explosive gases, or corrosive gases in the surrounding air.
  - Maintain a safety distance of 10 cm around the instrument while it is in operation.
- **Safety precautions:**
  - Do not unplug or plug in the power connector or toggle the power button with wet hands.
  - Do not unplug the power plug while the instrument is powered on.
  - Do not perform maintenance or cleaning while the instrument is powered on.
  - Do not install the instrument on an uneven, swaying, or vibrating work surface.