## Par de imanes cilíndricos, material AlNiCo

### Referencia QLL002

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
## Pair of cylindrical magnets, material AlNiCo

### Code QLL002

A magnet is a material or object that produces a magnetic field. This magnetic field is invisible but is responsible for the most notable property of a magnet: a force that pulls on other ferromagnetic materials, such as iron, steel, nickel, cobalt, etc. and attracts or repels other magnets.

A cylinder bar magnet is a type of bar magnet that is a cylindrical shape.


**Experiment 1:** To identify the magnetic and nonmagnetic substance.

#### Components Required:

- Magnet
- Various material (not included)

#### Procedure:

1. Test various objects of different materials to see if they are attracted by a magnet.
2. Test whether the two ends of the magnet (each identified by colour) act in the same way on these objects.

3. Repeat the experiment with the green end of the magnet and write the results in the table in the report.

**Table:**

<b>Object</b>	<b>The magnet is attracted by the red end</b>	<b>The magnet is attracted by the green end</b>
Glass Rod		
Iron Rod		
Aluminium		
PVC rod		

**Experiment 2:** To observe the magnetic field lines of cylindrical magnets.

**Components Required:**

- Cylindrical magnets
- Sheet of white paper (not included)
- Iron filing (not included)

**Procedure:**

1. Place the two magnets on a white sheet of paper so that same pole face each other.
2. Sprinkle some iron filing on the magnets.
3. Observe the magnetic field lines of magnets.