**SELLADORA PARA BOLSAS DE PAPEL PARA ESTERILIZAR
SEALING MACHINE FOR PAPER STERILISATION BAGS
MACHINE À SCELLER LES SACHETS DE STÉRILISATION EN PAPIER**

REF. - CODE - RÉF. - ZZC002

Mahita

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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TABLE OF CONTENTS

1. Technical specifications	4
2. Components	5
3. Operation	5
4. Maintenance	5

1. TECHNICAL SPECIFICATIONS

This equipment has been specially designed to ensure a firm and reliable seal on pouches and thermoplastic materials used in sterilisation processes. Its main function is to create an airtight seal that prevents the entry of contaminants and preserves the sterility of the contents. The sealer operates by the application of controlled heat, ensuring that the seal is uniform and durable, which is essential to maintain the integrity of the material until use. It is also designed to accommodate different bag types and material thicknesses, providing versatility and efficiency in the process of sealing packaging used in medical and laboratory environments.


Characteristics:

- Made of aluminium, with 4 rubber feet at the base
- Easy to operate, just press and release. Manual cutting by sliding blade.
- Reduced energy consumption because energy is only consumed when the sealing arm is pressed.
- Adjustable control with 8 levels to regulate the heating time in seconds according to the thickness of the material used
- Red indicator light during the warm-up time
- Includes a spare blade, 2 PTFE cover cloths and 2 heating elements.

Reference	ZZC002
Sealing width	3 mm
Max. sealing length	300 mm
Max. plastic thickness	0.4 mm
Dimensions	440x85x245* mm
Consumption	550 W
Weight	3.50 kg

*height with the arm released

2. COMPONENTS

1. Arm
2. Control knob
3. Plug

Warning: (1) This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. (2) Children should be supervised to ensure that they do not play with the appliance. (3) If the mains cord is damaged, it must be replaced by qualified personnel to avoid a hazard.

3. OPERATION

- Plug the cable into the mains socket.
- The equipment does not have a switch, it will only operate by pressing the arm.
- Set the control knob according to the thickness of the material to be sealed; use a higher number the thicker the material.
- Press the arm to start. Sealing is done when the red LED light is on. For better sealing, leave it pressed for one more second after the light turns off.
- Use a lower number on the control knob in case the seal burns. Increase the number if the seal is not airtight. If the bag sticks to the silicone rubber, you have not given enough cooling time.
- Once the time has been set, there is no need to readjust it as the equipment will remain constant and unchanged.
- Use the sliding blade to cut.
- The equipment will only consume power when the arm is pressed.

4. MAINTENANCE

- Unplug the power cord when servicing the equipment.
- Always keep the sealing platform clean. Leaving residues will reduce the life of the equipment.
- Never use moisture to clean the sealing surface.
- For optimum performance, use only original spare parts.
- A worn silicone rubber will influence the sealing, replace the rubber when it is worn or burnt.