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## **MAGNETIC FIELD CHAMBER**

**Article code: QLL007**

Magnetic field chamber demonstrates three-dimensional nature of magnetic field associated with bar magnets. It comprises a hollow chamber of clear acrylic with a tunnel at its center to hold cylindrical bar magnets. The chamber has high quality iron filings, which aligns itself as per the magnetic field geometry of a magnet. It also allows magnets of other shapes to be investigated by simple positioning their poles adjacent to any side of the chamber.

### **Procedure:**

1. Shake the magnetic field chamber and bubble from side to side to disperse the filings. Owing to the trapped pocket of air, the iron filings are evenly distributed within the vessel.
2. Insert the bar magnet through the center of the hollow axial bore and observe the changes.
3. Now without changing the position of the bar magnet, rotate the magnetic field chamber slowly till all the iron filings have attached themselves around the bar magnet.
4. Observe the alignment of iron filings.
5. Repeat the experiment with other magnets.

**Observations:** This demonstration reveals the following key characteristics of the magnetic field lines of a bar magnet:

1. Magnetic monopole does not exist, only dipoles and combinations of dipoles exist, which shows that the magnetic field lines begin from one pole and end at opposite pole. The magnetic field vectors are tangent to the lines at all points.
2. The density of magnetic field lines is proportional to the magnetic field strength in that region.
3. Magnetic field lines never cross themselves otherwise the field would have an ambiguous direction at the crossings.