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DESCRIPTION

Copper wire 16 SWG, mounted in U shape, vertically on the base. Maximum current is 8A.


INTRODUCTION

The electromagnetic field is a property of space caused by the motion of an electric charge. A stationary charge will produce only an electric field in the surrounding space. If the charge is moving, a magnetic field is also produced. An electric field can be produced also by a changing magnetic field. The mutual interaction of electric and magnetic fields produces an electromagnetic field, which is considered as having its own existence in space apart from the charges or currents (a stream of moving charges) with which it may be related. Under certain circumstances, this electromagnetic field can be described as a wave transporting electromagnetic energy.

MAGNETIC FIELD PATTERNS ASSOCIATED WITH A RECTANGULAR WIRE

The magnetic field lines around a straight conductor carrying current are concentric circles whose centers lie on the wire. The direction of magnetic field lines can be determined using right-hand thumb rule.

OTHER COMPONENTS REQUIRED

- Iron filing
- Power supply
- Patch cord

PROCEDURE

1. Connect the power supply with the vertical wire through the banana sockets.
2. Scatter iron filings on the stand of the vertical wire in such a way that only the near surface around the vertical wire is covered.
3. Switch on the power supply and gently tap the stand.
4. You will observe the ring-shaped magnetic lines of force around the wire.
5. As soon as you observe the magnetic lines of force, switch off the power supply.

In this experiment using vertical wire, you will observe ring shaped magnetic field lines around the vertical wire, which become weaker with increasing distance. According to Fleming's right-hand rule if the thumb of your right hand is pointing up in the direction of the current, then your fingers point in the direction of the magnetic field.