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## Magnetic field demonstrator, long solenoid

Article code: QLL010

### Description:

For the study and demonstration of magnetic field associated with a current carrying long solenoid. Mounted on clear, transparent acrylic permits its use on OHP for classroom demonstration. Accompanied with a pair of 4mm socket terminals for electrical supply. The magnetic effects can be visualized through the use of iron filings or plotting compasses. DC power supply of 0-25volt/0-20 ampere current is required.

**Long solenoid:** A solenoid coil about 52×125mm (diameter × length) of 16 SWG enameled copper wire. Each turn is sufficiently apart and well insulated from the adjacent ones.

### Introduction:


**Electromagnetic field**, a property of space caused by the motion of an electric charge. A stationary charge will produce only an electric field in the surrounding space. If the charge is moving, a magnetic field is also produced. An electric field can be produced also by a changing magnetic field. The mutual interaction of electric and magnetic fields produces an electromagnetic field, which is considered as having its own existence in space apart from the charges or currents (a stream of moving charges) with which it may be related. Under certain circumstances, this electromagnetic field can be described as a wave transporting electromagnetic energy.

### Magnetic field patterns associated with Solenoid:

Solenoid consists of a length of insulated wire coiled into a cylinder shape.

- Current in solenoid produces a stronger magnetic field inside the solenoid than outside. The field lines in this region are parallel and closely spaced showing the field is highly uniform in strength and direction.
- Field lines outside the solenoid are similar to that of a bar magnet, and it behaves in a similar way – as if it had a north pole at one end and South Pole at the other end. Strength of the field diminishes with distance from the solenoid.
- Strength of the magnetic field can be increased by:

- a) Increasing the current in the coil
  - b) Increasing the number of coils in the solenoid
  - c) Using a soft iron core within the solenoid
- Reversing the direction of the current reverses the direction of the magnetic field.

#### Other components required:

- Iron filing
- Compasses
- Power supply
- Patch cords

#### Working:

1. Connect the power supply with the solenoid through the banana sockets.
2. Scatter iron filings on the stand of the solenoid in such a way that only the approximate surface around the solenoid is covered.
3. Also place the magnetic compasses near the solenoid.
4. Switch on the power supply and gently tap the stand.
5. You will observe the ring-shaped magnetic lines of force around the solenoid.
6. As soon as you observe the magnetic lines of force, switch off the power supply.
7. In this experiment using solenoid, several current carrying conductor loops are added together.

According to Fleming's right hand rule if the thumb of your right hand is pointing up in the direction of the current then, your fingers point in the direction of the magnetic field.