**AGITADOR VORTEX MODELO 740**  
**VORTEX AGITATOR MODEL 740**  
**AGITATEUR VORTEX MODÈLE 740**

Ref. | Code | Réf. LGB009

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*This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.*

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**TABLE OF CONTENTS**

INTRODUCTION .....	4
APPLICATION .....	4
TECHNICAL SPECIFICATIONS .....	5
OPERATION .....	5
<i>Instrument Operating Environment</i> .....	5
<i>Safety Precautions</i> .....	5

**INTRODUCTION**

The Nahita-Blue Model 740 vortex agitator incorporates international manufacturing technology and advanced design, combining functionality with a modern and elegant aesthetic. Its compact and robust structure makes it highly versatile, while its easy operation ensures a comfortable and efficient user experience.

It is specially designed for the effective mixing of a wide variety of samples that require thorough homogenization.

**APPLICATION**

It is widely used as a tool in research and experiments at universities, in the environmental protection sector, public health, disease prevention, the chemical industry, drinking water treatment plants, and in a variety of other applications.

## TECHNICAL SPECIFICATIONS

<b>Code</b>	<b>LGB009</b>
Model	Model 740
Control	Constant speed
Base material	Zinc alloy
Motor type	Brushless motor
Orbital diameter	4.5 mm
Maximum agitation volume	50 ml
Speed	3000 rpm
Operating conditions	5-40 °C / 80% HR
Protection level	IP42
Power	12 W
Power supply	100-240 V 50/60 Hz
Dimensions	93x93x72 mm
Weight	0.40 Kg

## OPERATION

- Place the sample to be mixed at the center top and lightly press to initiate the agitation process.

### Instrument Operating Environment

- This device is designed for indoor use only
- Altitude: ≤2000 m
- Operating temperature range: +5 °C to +40 °C
- Humidity range: ≤80%
- Place the device in an environment free from vibrations or air currents that may affect its performance
- Do not use the device in environments with conductive dust, explosive gases, or corrosive vapors in the surrounding air
- Maintain a safety clearance of 10 cm around the instrument while in operation

### Safety Precautions

- Do not plug or unplug the power cord, or operate the power switch, with wet hands
- Do not unplug the power cord while the instrument is powered on
- Do not perform maintenance or cleaning while the instrument is on
- Do not install the instrument on an uneven, unstable, or vibrating surface