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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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
2. Operation

1. Connect the pump to the power source, turn on power switch.
2. Connect the pump and the vacuum equipment with a hose.

3. Use criteria

1. Mounting the pump on a stable, rigid operating surface.
2. Remove the negative pressure before restart if stop the pump in process.
3. Make sure there is no suck the solvent or impurity into the pump.
4. Stop the pump in time when the process is over. The pump will switch to guard mode when the motor is too hot, restart when the temperature get down.

4. Troubleshooting

1. If pump fails to start, check the pressure gauge whether the equipment has a negative pressure, remove the negative pressure and restart.
2. If the motor fails to start, check the power and fuse.
3. Maintenance. If the solvent is sucked into the pump or the pressure cuts down, repair it as follow:
 - Remove the six screws in the bottom of the pump, pull the handle and remove the shell. Do not break the wire in the pump.
 - Remove the four stainless steel screws on th head of pump, remove the cover of head, valve plate, diaphragm one by one, clean them carefully and them fix them.
 - Tightening the four stainless steel screws on the head to deal with the pressure cut down.

 Note: The period of service time is one year since be sold. The parts and service are free except break by wrong operation.

Environmental protection instructions

Do not dispose of this equipment with normal household waste at the end of its life cycle; take it to a collection point for the recycling of electrical and electronic equipment. It does not contain elements that are hazardous or toxic to humans but improper disposal will harm the environment. The materials are recyclable as indicated by the marking. By recycling materials or other forms of reuse of old equipment, you are making an important contribution to the protection of the environment. Please contact your local authority for advice on the collection points.

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
CONTENT INDEX

1 Principle and characteristics about the pump	4
1.1. Characteristic	4
2 Operation	5
3 Use criteria	5
4 Troubleshooting	5

1. Principle and Characteristic about the Pump

Oilless diaphragm vacuum pump is the important equipment of enhance the filtration efficiency effect. Usually, make use of oilless diaphragm vacuum pump in order to prevent the medium be polluted. The pressure of the pump in round numbers 0.08Mpa. The pump's flux is a key to enhance filtration speed, especially to a large number of medium or the medium is too thick. The service life of the pump depends on the right operation, a reasonable motor power and the format of diaphragm.

1.1. Characteristic

HP series are designed for the Solvent Filtration. It is small, practical, and efficient. It is also easy to operate with low noise.