**INCUBADORAS REFRIGERADAS BIOQUÍMICAS
BIOCHEMICAL REFRIGERATED INCUBATORS
INCUBATEURS RÉFRIGÉRÉS BIOCHIMIQUES**

Ref. JBJ004 - JBJ005

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6. MAINTENANCE

- The water accumulated in the plastic tank must be discharged regularly.
- This incubator should be used indoors, keeping ventilation or air conditioning, otherwise, the space sealing will cause the temperature in the room to rise because of the heat dissipation of the incubator, it will make the working environment worse, and the high ambient temperature will affect the normal operation of the device.
- The distribution equipment of the user shall meet the working needs of incubator. In order to ensure the safety of the equipment and test, please install the external leakage protection device and supply the power supply according to the requirements of the equipment nameplate. The power supply line in the workplace must meet the demand of the corresponding power, and the power supply must have reliable grounding.
- Clean the condenser every two months, and avoid the decay of the cold efficiency of the compression mechanism.
- When you stop working, turn off the incubator. When you do not use the device for a long time, unplug the power supply and cover the equipment with soft cloth or plastic paper to prevent dust from entering.
- If the power cord is broken, it must be replaced. The replacement must be replaced with the same type and specification of power cord.
- The power plug of the incubator should be unplugged from the power outlet immediately and the supplier shall be contacted or the trained maintenance personnel shall be asked to deal with it in the following cases:
 - The liquids drop into the inside of the device.
 - The equipment is drenched or watered by rain.
 - Falling of the device, strong vibration or damage to the shell.
 - The functioning of the incubator has changed obviously.
 - Abnormal sound or smell occurs when the equipment is not working properly.
- The user should not try to repair the device, which will disqualify you from warranty, if it is in validity period, and may be subjected to electric shocks. For repairs, please contact your distributor.

Instructions on environment protection

At the end of its life cycle, please do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans, but a non adequate disposal would damage the environment. The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment. Please inquire at the community administration for the authorized disposal location.

INDEX OF LANGUAGES

Spanish	2-10
English	11-18
French	19-26

INDEX OF CONTENTS

1. Safety warnings	11
2. Product overview	12
3. Technical parameters.....	13
4. Installation instructions.....	13
5. Control panel.....	14
6. Maintenance	18

1. SAFETY WARNINGS

PROHIBITIONS

- Do not place volatile, flammable or explosive materials near this equipment, otherwise it may cause explosion or fire.
- Do not place the equipment in rain, humidity or water drenching places, otherwise it may lead to electric leakage, short circuit or electric shock.
- Do not insert metal objects, such as iron wires into any orifice and gap of the device, in order to avoid electric shock or injury.
- Non-professional persons are not allowed to disassemble, repair or modify this equipment, otherwise it may cause fire or electric shock accident due to improper operation.
- Do not damage the power plug or power wire. If the plug is loose, the power cord must be replaced.

OBSERVING ITEMS

- This equipment should be installed on solid ground. If the ground is not solid or the installation place is not appropriate, the person may be injured due to the device being tipped over.
- Please use the power supply indicated on the nameplate of this equipment. The equipment must be reliably grounding in case of accidental leakage of electricity or fire.
- Before any repair or maintenance of this equipment, it is important to disconnect the power supply to prevent electric shock or personal injury.
- Be sure to wear gloves for equipment repair or maintenance to prevent damage caused by touching edges or sharp corners.
- Do not damage the power wire or use an unspecified power wire, otherwise it may lead to electric shock or fire.
- If an abnormal operation is found, power off immediately and stop the operation. It is possible to cause electric shock or fire in the abnormal state.

ATTENTION ITEMS

- The equipment should be placed on a horizontal surface.
 - Use an independent power socket with ground wire and insert the plug tightly.
 - Unplug the power cord before moving the incubator to a new location.
 - Internal parameters setting should be responsible by the specific management personnel, in order to prevent the controller program dysfunction.
 - Equipment installation location should have good ventilation. Should be placed away from the wall;
 - the distance from objects must be greater than 30 cm.
- The surfaces shall be not contacted with the volatile chemicals such as gasoline, thinner, etc.
- Keep the incubator inside and outside clean.

2. PRODUCT OVERVIEW

The biochemical refrigerated incubator is a device widely used in environmental protection departments, health and epidemic prevention departments, agriculture research institutes, etc. Can be used for water analysis, BOD determination (BOD stands for Biochemical Oxygen Demand), microbial culture, breeding experiments, etc.

Main characteristics

- The shell is made of cold rolled steel plate, and the working chamber and shelves are made of mirror stainless steel, convenient for cleaning and disinfection.
- Casters are provided at the bottom for easy movement; the two front casters have brakes.
- The air-inlet grille can be opened directly to facilitate the cleaning of the condenser (please clean the condenser every two months so as not to affect the refrigeration performance of the compressor).
- Control panel with color LCD screen and 4 buttons.
- Over temperature alarm. Additional temperature limiter.
- High quality compressor ensures long running time.
- Double evaporator switch ensures no frost operation.
- Special design ensures temperature uniformity within 1°C.

6. Compressor working mode and defrosting related parameters

■ Press the SET key for about 3 seconds, and the password prompt “Lc” appears. The password is modified to “123” by adding and reducing keys. Click the setting key to enter the compressor setting state, and then click the setting key for 3 seconds, you can exit this state, the parameters are automatically saved. If there is not any bottom within 30 seconds, you can automatically exit this state, and the current parameter is not saved.

Parameter indicator	Parameter function description	(Setting range) factory set value
Lc-	You can access to and modify parameter values when Lc=123	
S-H	When “S-H> temperature set value”, the compressor works continuously; otherwise, the compressor works intermittently. Note: when “S-H=50.0”, the compressor does not work	(-20.0~50.0)
Ft-	After defrosting, the fan delays Ft- to turn on again	(0~99s) 30s
dt1	Defrosting interval 1 when SP ≤15°C	(0~250 horas) 6 horas
HS1	Defrosting output 1	(0~250s) 60s
dt2	When 15 °C<SP≤30 °C defrosting interval 2	(0~250 horas) 12 horas
HS2	Defrosting output 2	(0~250s) 60s
HA-	The compressor is turned on and off automatically; 1 is to automatically judge the switch compressor according to the ambient temperature, and 0 is to control the compressor to open and close according to uP and dn.	(0~1)1
cH-	When the defrosting and evaporator switch the selection; 0 is defrosting, 1 is evaporation switching.	(0~1)0
AP-	Compressor protection; if set value ≥ AP, the compressor will not open.	(0~AP) 50.0
bP-	When intermittent operation and the set value ≥ bP + ambient temperature, the compressor will not open.	(0~AP) 6.0

7. Notes

- Working mode of compressor.
 - After the electric delay of 30 seconds on the controller, the compressor can be automatically judged in what way it works.
 - When the temperature setting value changes beyond 3.0 °C, the controller automatically determines how the compressor works. If the set value is not greater than 3.0 °C, the compressor will remain in its original mode.
 - When the set values of “S-H” and “HA” are changed, the “S-H” parameter must be re-powered to the controller.
- Warning function
When the controller is powered on for the first time or the temperature setting value is reduced, even if “PV>=SP+AL”, the alarm light is not on and the buzzer does not ring.
- Temperature failure
If the temperature measurement shows “oooo”, indicating a temperature failure, please check the sensor first. If the sensor has no fault, the controller fails.

4. Contact address

■ Press the SET key for about 3 seconds, and the password prompt “Lc” appears. The password is modified to “9” by adding and reducing keys. Click the setting key to enter the address setting state, and then click the setting key for 3 seconds, you can exit this state, the parameters are automatically saved. If there is not any bottom within 30 seconds, you can automatically exit this state, and the current parameter is not saved.

Parameter indicator	Name	Parameter function description	(Setting range) factory set value
Lc-	Password	You can access to and modify parameter values when Lc=9	
Ad-	Contact address	Local communication address	(1~32)1
rHd-	Ambient temperature correction		(-20.0~20.0 °C) -2.0 °C
ET	Timing function	When ET=0, no timing function; when ET=1, as soon as the incubator is turned on, the timer begins; when ET=2, when the temperature reaches the set point, the timer begins	(1~2) 2
En	End of operation thermostat	When En=0, when the operation is over, stop heating When En=1, when the operation is over, continue to maintain a constant temperature	(0~1) 0
rL-	Range value	The temperature setting lower limit	(-15~50 °C) 0 °C
rH-	Range value	The temperature setting upper limit	(0~95 °C) 70 °C
P-t	Print interval	When P-t=0, it will not print.	(0~9999) 0
Lt-	Output power	Percentage of output power	(0~100%)100
co-	Turn-off deviation	When the measured value>set value+co, the heating output is turned off	(0~20.0 °C)1.0

5. Environment temperature

■ Press the SET key for about 3 seconds, and the password prompt “Lc” appears. The password is modified to “18” by adding and reducing keys. Click the setting key to enter the environment temperature setting state, and then click the setting key for 3 seconds, you can exit this state, the parameters are automatically saved. If there is not any bottom within 30 seconds, you can automatically exit this state, and the current parameter is not saved.

Indicador parámetro	Name	Parameter function description	(Setting range) factory set value
Lc-	Password	The environment temperature is shown in “Lc=18”	
Hd	Current ambient temperature		
HE-	Ambient temperature measurement enabled	0:Hd= HSET 1:Hd=measured value	(0,1) 1
HSET	Ambient temperature setting value		25.0

3. TECHNICAL PARAMETERS

Code	JBJ004	JBJ005	
Temperature	Temperature range	0~65	
	Temperature fluctuation (°C)	±0.5	
	Temperature resolution (°C)	0,1	
	Temperature accuracy (°C)	1°C (a 37°C)	
	Working ambient temperature (°C)	+5~35 °C	
Structure	Rated power (W)	600	800
	Compressor	Air-cooled hermetic compressor	
	Refrigerant	R134a	
	Defrosting	Double evaporator switch	
	External power supply	One five-hole socket (on the left side of the incubator)	
	Testing hole	Inner diameter of 52 mm (on the left side of the incubator)	
	Insulation material	Sponge	
	Controller	Control panel	Color LCD screen with 4 buttons
Temperature control		Single stage timing control	
Sensor		PT100	
Additional function		Fixed value mechanical temperature limiter	
Other specifications	Internal size (width*deep*high, mm)	420*355*550	500*405*750
	External size (width*deep*high, mm)	588*670*1210	668*720*1410
	Capacity (L)	70	150
	Max. shelf load (Kg)	15	
	Net weight, (kg)	85	120
	Power supply	AC220V, 50Hz	
	Shelves included	2	
	Optional capabilities	RS485, USB port	

4. INSTALLATION INSTRUCTIONS

- After opening the package, check that the device is intact, and the incubator is placed on the flat floor of the room.
- The equipment should be placed in the place of air circulation, no vibration, no flammable, explosive gas or high dust. There is no conductive dust, corrosive gas and no strong magnetic field interference around the incubator. Should avoid direct sunlight.
- In order to ensure the safety of equipment and test, please install external leakage protection device

and supply the power supply according to the requirements of equipment nameplate.


- The blue PU flexible tube must be inserted into the drain hole of condensate water, and the condensate water is connected with the water receiver.
- In the boot step, plug in the power supply to turn on the main power switch. At this point, the device is lit, and the corresponding temperature can be set as needed to start running.
- In the turn off step, the incubator is out of power and stops working, you can unplug the plug.

5. CONTROL PANEL

I. Panel indications

Display indicators

- | | | | |
|--|----------------------------------|--|-----------------------------|
| | Temperature setting icon | | Alarm icon |
| | Time setting icon | | Refrigeration icon |
| | Internal parameters setting icon | | Heating icon |
| | Running status icon | | User parameter setting icon |
| | | | Fan icon |
| | | | Stop running status icon |

2. Operation method

■ Click on the setting key, it will enter the temperature setting state. The required value is modified by increasing and reducing keys. Press the setting key, it will enter the time setting state; the time hour area is flashed, then press the setting key, the time minute area flashes, which can modify the required value by increasing and reducing keys, and press the setting key again to save and exit the setting state. When the time is set to "0", it indicates that there are no timing function. When the setting time is not "0", the measuring temperature reaches the setting temperature. The timer begins to clock, the time is up, and the run is over. "Stop" character lights up, the buzzer calls 30 seconds (any key can be used to

silence). If you press the reducing key for 4 seconds, the program starts running again.

- The chamber lighting may be turned on or off at the customer's request, with a corresponding light on or off.
- Click on the "Add" key to switch the LCD backlight in an unset state.

3. Temperature internal parameters

■ Press the SET key for about 3 seconds, and the password prompt "Lc" appears. Through adding and reducing keys, the password is changed to "3". You can click the setting key to enter the temperature internal parameter state. Click on the setting key to modify the parameters, and then press the key for 3 seconds to exit this state and save the parameters automatically. If the SET key is not pressed within 30 seconds, this state is exited automatically and the current parameter is not saved.

Note: the internal parameters have been adjusted before leaving the factory. Please do not adjust without special needs.

Parameter indicator	Name	Parameter function description	(Setting range) factory set value
Lc-	Password	You can access to and modify parameter values when Lc=3	
AL-	Over temperature deviation alarm	When "PV>=SP+AL", the buzzer rings and has a refrigeration output	(0.0~20.0) 5.0 °C
cT-	Compressor start-up delay	The compressor delays to protect time, twice start time ≥ cT seconds	(0.0~600 segundos) 180
T-	Control cycle	Heating control cycle	(1 - 60 segundos) 5
P1-	Proportional band	Normally open time scale action regulation	(1.0~valor rango) 20
I1-	Integration time	Normally open integral action regulation	(0-1000 segundos) 200
d1-	Rate time	Normally open differential action regulation	(0-1000 segundos) 200
P2-	Proportional band	Disconnected time scale action regulation	(1.0~valor rango) 35
I2-	Integration time	Disconnected integral action regulation	(0-1000 segundos) 500
TOF-	Deviation temperature	Disconnect working deviation temperature	(0~10.0) 3.0
uP1-	Compressor start threshold 1	The compressor is in manual mode and uP/dn is controlled by uP1/ dn1 when the set value < ambient temperature - deviation temperature. uP/ dn is controlled by uP3/dn3 when environment temperature - deviation temperature <= set value < environment temperature+ deviation temperature. uP/dn is controlled by uP3/dn3 when set value >= environment temperature+ deviation temperature. When "PV>=SP+uP" and the compressor delays the time, the compressor is opened.	(-4.9~5.0) 0.1
dn1-	Compressor close threshold 1		(-5.0~uP-0.1) 0
uP2-	Compressor start threshold 2		(-4.9~5.0) 0.3
dn2-	Compressor close threshold 2	The compressor is in manual mode, and when "PV<=SP+dn", the compressor is closed.	(-5.0~uP-0.1) 0.2
uP3-	Compressor start threshold 3		(-4.9~5.0) 0.4
dn3-	Compressor close threshold 3		(-5.0~uP-0.1) 0.3
Pb-	Zero setting	Zero error correction of sensor Pb=actual temperature value - instrument measured value	(-9.9~9.9) 0.0
PK-	Full adjustment	Sensor full scale error correction PK=1000*(actual temperature value - instrument measured value)/instrument measured value	(-999~999) 0

Note:

1. In order to achieve the best effect of the control, P can be modified appropriately (the modification range is generally between "10.0 and 20.0"), but the factory parameters can usually be used directly without modification.