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- Keep the oil clean. If it becomes dirty, cloudy, or has water or other volatile substances that affect the limit vacuum, the oil should be changed. Before changing the oil, start the pump and let it run for about 30 minutes to dilute the oil, then stop the pump and let the oil flow out of the oil outlet. Then run it for 1-2 minutes with the gas inlet open. During this time, add a small amount of clean oil through the oil inlet to replace any remaining oil in the pump. Do it repeatedly. After making sure the pump is clean, fit the oil drain plug.
- If it is not going to be used for a long time, cover the grease plug and the air plug and store it in a dry place. In the meantime, protection against moisture and rust must be taken into account.
- If it is necessary to disassemble and repair, make sure that an experienced or professional technician does it.

## 7. Warranty

- Within 14 days from the delivery date if manufacturing defects are detected, you can return the item to us for replacement, attaching the purchase invoice. But please make sure there is no physical damage to the equipment and packaging before you replace the item.
- If after 14 days from the delivery date the item has a defect, you can call the local distributor for your claim.
- The vacuum pump has a 1-year warranty for non-wear parts and the motor, except for physical damage due to misuse of the same

### Environmental protection instructions

Do not dispose of this equipment with normal household waste at the end of its life cycle; take it to a collection point for the recycling of electrical and electronic equipment. It does not contain elements that are hazardous or toxic to humans but improper disposal will harm the environment. The materials are recyclable as indicated by the marking. By recycling materials or other forms of reuse of old equipment, you are making an important contribution to the protection of the environment. Please contact your local authority for advice on the collection points.

## INDEX OF LANGUAGES

Spanish .....	2-5
English .....	6-9
French .....	10-13

## CONTENT INDEX

1 Technical specifications .....	6	6 Maintenance .....	8
2 User's manual .....	7	7 Warranty .....	9
3 Cautions .....	7		
4 Installation .....	7		
5 Troubleshooting guide .....	8		

## 1. Technical specifications

Code	PXB012	
Voltage	110V / 60Hz	220V / 50 Hz
Free air displacement	5 CFM	4,5 CFM
Full vacuum	3x10 <sup>-1</sup> Pa	
	0,05 Mbar	
Rotation speed	3500	2800
Power	1/2HP	
Oil capacity	330 mL	
Dimensions	290x115x220 mm	
Weight	9 Kg	

## 2. User's manual

1. Check the oil level before operation, making sure that it is not below the level line. Be sure to fill the oil with the correct level and type of oil before starting (use high-speed vacuum oil type "HFV 32")
2. Remove the air fill cap and connect the pump canister. The tube must be short and sealed to prevent any leakage.
3. Remove the exhaust plug, plug in and turn on for operation and use.
4. After using the machine, unplug it, remove the connection pipes and plug the exhaust and oil plug.

## 3. Cautions

- Do not pump inflammable, explosive and poisonous gases.
- Do not pump gases which can corrode metals and can exert chemical changes with the pump oil.
- Do not pump gases containing pellet dust.
- The temperature of the pumped gas should not be over 80°C and the environment temperature should not be lower than 5°C – 6°C.
- Do not use it as a compression or a convey pump.
- Do not use or operate without oil.
- Do not pull on the wire when unplugging from the outlet.
- Do not put heavy weights on the wire to avoid being extruded..
- Do not use when the plug or outlet if it is damaged
- Do not plug or unplug with your wet hands.
- Unplug or switch off the machine when there is any kind of leak of coal gas.

## 4. Installation

1. When using the pump, it must be in a horizontal position and place it in a dry and clean place with adequate ventilation. The distance from the periphery must be more than 2 cm and the distance between the front and rear must be at least 5 cm. When installing the equipment, make sure that a enough amount of air can enter through the housing side of the paddle.
2. Installation in a complete set:
3. Remove the rubber cap from the bottom base and fix with the screw. If you do not have special equipment for installation, we invite you to be assisted by our Technical Service for this task.
4. If the gases from the pump are harmful to health or have a bad effect on the environment, you can extend the pipes to let the exhaust mouth go outside, however, be sure to practice a logical method that protects the environment.

## 5. Troubleshooting guide

FORM OF BREAKDOWN	CAUSES	REMOVING WAYS
LOW DEGREE OF VACUUM	1. Lack of oil	Fill the oil at the above level line
	2. The pumping oil is emulsified or dirty	Change the oil
	3. The oil-entering mouth is blocked up short in oil	Clean the oil-entering mouth, wash the filter net
	4. The connection of drain pipes	Check the connection of the pipes
	5. The pump is not suitable	Verify the proper use of this pump
OIL DRAIN	1. The oil-selter is damaged	Change it
	2. The tank is loose or worn out	Screw it and the O shaped ring
OIL SPRAY	1. Too much oil	Put the oil at the level line only
	2. The pressure of the gas-entering mouth is too high or it has pumped much	Choose a bigger pump
DIFFICULT IN STARTING	1. The oil temperature is too low	Starts the machinery several times or heat the oil
	2. There is some breakdown in the electrical machinery	Check and fix it
	3. Some foreign matter is in pump	Check to remove

## 6. Maintenance

- Keep the pump clean and avoid the entry of foreign bodies.
- Maintain the oil level. Don't let it run without oil.