MICROSCOPIOS SERIE P
MICROSCOPES SERIES P
MICROSCOPES SÉRIE P

Zuzi

Ref. HBC023 - HBC024

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If the field of vision is cut:

Check the nosepiece is properly fit. To do that, turn it slightly back and forth until it fits.

If there is dust or dirt in the field of vision:

Check there is dust on the precondenser's lens, upper lens of the condenser or the eyepiece, as well as the cleanliness of the sample. Once you find out where the dirt is you should clean it as explained before.

If the sample gets fuzzy while observing:

Check the tension adjusting ring, as it may be too weak and provoke an involuntary fall of the stage.

If an area in the field of vision is out of focus:

Check the objective is properly located on the luminous beam's path.

If not you should turn it until it is appropriate position.

Check the sample is properly located on the stage, making sure it is held by the stage clips. Then, if it does not focus you should revise the pre-focusing control.

If the frontal lens of the objective touches the prepared slide when focusing:

Check the prepared slide is not the other way round (being the cover glass over the slide) and place it properly in case it is appropriate to do so. Check also the slide's thickness, which has to have a standard thickness of 0,17 mm.

8. RECOMMENDATIONS

When using the oil immersion objective you should notice the following:

- In order to use all the numerical aperture of the oil immersion objective:
 - Pre-focusing: Focus the sample with the objective with less magnification (4x).
 - Put a drop of oil immersion on the cover glass.
 - Rotate the nosepiece, place the immersion objective and arrange the focusing with the fine adjustment knob.
- You should avoid air bubbles in the oil so as not to spoil the resulting image. To do so you should move the pre-focusing control up and down to make the oil layer uniform.
- Turn the nosepiece and then locate the oil immersion objective; you can adjust focusing with the fine knob.
- We highly recommend drying the frontal lens after each use. You should never use a duster. Instead of it you should use specific optics cloth and liquids that are sold to this purpose (as the ones for cleaning sunglasses or binocular lenses). You should not submerge the objective in any liquid, as this could provoke the lens' detachment. Depending on the observer preferences the tension of the coarse adjustment knobs can be modified, since it is provided with a tension adjustment ring that is located next to the right coarse knob of the microscope. However, you should avoid loosening it excessively, as this could produce an involuntary fall of the stage and so provoke accuracy deterioration.

9. ACCESSORIES

- Power cord
- Dust cover
- Blue filter (32mm)
- Spare fuse 2A (5x20mm)
- User manual
- Dipping oil immersion if the microscope incorporates the 100x lens.

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INDEX OF CONTENTS

1. Application	11	5. Maintenance	15
2. Description	11	6. Causes of a defective image.....	15
3. Technical specifications.....	13	7. Troubleshooting.....	15
4. Installation.....	13	8. Recommendations.....	16
4.1 Instrument installation	13	9. Accessories.....	16
4.2 Using the instrument	13		

Thank you for choosing this equipment. We sincerely wish that you enjoy the microscope and we highly recommend looking after this equipment according to what is stated in this manual. Read the point of Maintenance where we expose the security notes.

Please, bear in mind the following:

- This manual is inseparable from the microscope so it should be available for all the users of this equipment.
- It should be handled carefully avoiding sudden movements, knocks, free fall heavy / Sharp object on it.
- Avoid spilling liquids inside the equipment.
- To prevent fire or electrical discharges avoid dry or dusty environments. In case it may happen, unplug the equipment immediately.
- If you have any doubt (installation, setting up, operation), please contact your distributor.
- Never dismantle the different pieces of the microscope to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the microscope.
- Operations made by non-qualified staff will automatically produce a loss of the microscope warranty.
- Neither fuses nor accessories are covered by the warranty. The warranty neither covers piece's deterioration due to the course of time.
- Overhaul is not covered by the microscope warranty.
- Make sure you keep the invoice, either for having the right to claim or asking for warranty coverage.
- Manufacturer reserves the right to modify or improve the manual or equipment.

5. MAINTENANCE

- The microscope must be installed in a dry, clean and free of acid and alkaline vapors space; avoid the direct exposure to direct sunlight.
- Environmental conditions: the temperature must be between 0°C- 40°C and the relative humidity at maximum of 85%.
- The microscope has been adjusted and inspected strictly before leaving factory; handle it carefully avoiding hits and knocks.
- If there is dust in the lenses blow it with a rubber blower and then clean softly the lens with a soft brush; clean carefully the rest of oil or finger print from the Surface of the lens with a soft cloth moistened with a few organic solvent (mixture of ether and alcohol 7:3).
- Do not wipe the lens surface very regularly, or else the lens will be scraped, reducing the quality of the transmission and imaging. Please keep the rest of the instrument clean.
- Keep the mechanical parts clean and dry.
- Unplug the power cable and disconnect the microscope from the power supply in case of not being used for a long time; adjust the light intensity to the minimum and cover the microscope with the dust cover.
- To replace the fuse, open the fuse holder with a flat screwdriver. Remove the fuse and replace it.

6. CAUSES OF A DEFECTIVE IMAGE

In case of a defective image you should check the following:

- Illumination is well done and luminous intensity is neither excessive nor too weak.
- Eyepieces are well fit and objectives well screwed.
- When cleaning the entire optical system you should proceed as follows:
 - Make the eyepieces turn checking the little specks are moving; if so you should clean them. Make the head turn. You should never disassemble the head, but you can clean it with delicacy by blowing the accessible surfaces of prisms with a plastic bulb.
 - The objective and, in case the parasite images also turn, clean it with the help of a dry brush to remove the dust. Watch the front part of the microscope with a magnifying glass or an inverted eyepiece.
- Immersion oil is enough and does not contain air bubbles or impurities.
- Thickness of slides, cover glasses and means of assembly is not too much as to avoid focusing with medium or high magnification (there are standard sizes for both the slides and cover glasses). Slides and cover glasses should be clean. Check that the slide is properly located and there are not two cover glasses superimposed.

7. TROUBLESHOOTING

Before sending the microscope back to Technic Technical Assistance Department you should check the following:

If the lamp does not switch on:

Check it is properly connected to the right current intake.
Check both fuse and lamp are in good conditions. If not you should replace the fuse.

Insert the plug into the socket at the back of the microscope. Connect the other end to the electrical power supply. If the batteries are charged the microscope can be used without connecting to the electrical power supply.

Note: in case of using the microscope connected to the power supply, it should be properly earthed. Make sure the input voltage is the appropriate for this equipment (100-240V; 50-60Hz).

4.2 Using the instrument

- Turn the power source of the microscope on by pressing the switch at the back of the microscope. The light intensity can be regulated by means of the knob placed at the left side of the base of the microscope. It is recommended to be used at 70% of intensity.
- Place the simple (microscope slide) and hold it with the clamp; to open the clamp use the metal lever at the left side.
- The amplitude of the light beam that reaches the simple can be adjusted by changing the diaphragm aperture. The maximum resolution is obtained by adjusting properly this aperture. Every time the objective is changed and to get the best resolution it is advisable to remove the eyepiece to see the size of the diaphragm aperture through the eyepiece tube. The diaphragm aperture must be adjusted until it is a little smaller than the objective aperture.

Note: The diaphragm aperture does not change light intensity; for this use the knob placed at the left side of the base of microscope.


- If needed put the filter in the filter holder.
- Adjust the interpupillary distance (distance between eyepieces) before starting observation. Moreover, it will be necessary firstly to focus image with one of the eyepieces and then adjust the diopter correction of the other eyepiece until the image is also focused with this eyepiece.
- Start observation always with low magnification objectives 4x and 10x. To avoid the objective touching the sample, turn the coarse focusing knob to lift the mechanical stage until it reaches a few millimeters from the objective. Then turn slowly the fine focusing knob to focus the image.
- Once the image is well focused change to the next higher magnification objective and adjust focus by turning only fine focusing knob. To change the objective turn the nosepiece until hearing a "click" sound and make sure it is well centered in the light beam. Always try to change the objective by holding the nosepiece not the objective itself; take also into account the height of each objective since the 100x objective is longer and it could touch the sample.
- When using 100x objective, move the condenser to the highest position and use a drop of immersion oil; avoid the formation of bubbles that could affect observation of sample. The 100x objective must be wiped off after each use with a lens paper to remove all the remains of immersion oil.
- When turning the coarse and fine focusing knobs to move the stage, if it is too hard or too loose it will be necessary to adjust the tension adjustment ring. The tension adjustment ring is placed in the right coarse focusing knob.
- Move the sample on the stage by using the translation knobs to study the different structures of the specimen.

1. APPLICATION

Zuzi microscopes P Series have been designed and developed to be used in different fields of applications as biology, bacteriology, histology or pathology. They can also be used in education and University research.

2. DESCRIPTION

The different parts of the microscope are shown on the below image (Figure 1)

Depending on the model the head can be binocular (HBC023) or trinocular (HBC024). The Siedentopf head is 30° inclined and 360° rotatable with a lateral screw to fix its position. The interpupillary distance can be easily adjusted (55-75mm) and they have diopter correction in both eyepieces (Figure 2). The trinocular HBC024 is provided with a height adjustable third tube in which a USB videocamera can be placed.

Figure 2

The stand makes it easy to carry and avoid damages when trying to hold it by other components. The microscope is provided with a quadruple inverted nosepiece to make changing the sample easier and more comfortable.

The stage is a double layer mechanical stage with translation coaxial knob and presents a double clamp to hold up to two slides. The microscope is provided with a safety screw to fix the maximum upper position of the stage and avoid the stage to hit and damage the objectives (Figure 3).

Figure 3

The coaxial fine and focusing knobs are at an ergonomic position and are equipped with a tension adjustment ring (Figure 4).

Figure 4

The microscope is provided with rechargeable batteries (Figure 5) to be used without connecting to the power supply. The batteries are charged when the microscope is connected to the electrical power supply.

ATTENTION!:: DO NOT USE THE MICROSCOPE WITH ALKALINE BATTERIES TO AVOID SECURITY PROBLEMS.

Figure 5

3. TECHNICAL SPECIFICATIONS

CODE	HBC023	HBC024
HEAD	Binocular	Trinocular
	Tipo Siedentopf; inclinado 30° y giratorio 360°	
	Distancia interpupilar regulable (55-75mm); ajuste de dioptrías en ambos oculares	
EYEPIECES	WF10x/18 mm	
NOSEPIECE	Quadruple inverted	
OBJECTIVES	Plan achromatic objectives: 4X, 10X, 40X (R), 100X (R) (I).	
STAGE	Double layer mechanical stage, Size 132x 142 mm Movement 75x50mm	
CONDENSER	Abbe NA 1.25 with iris diaphragm and filter	
FOCUS	Coaxial coarse and fine focusing knobs; focus scale 0.002mm	
LIGHT	LED 3W/4V, intensidad regulable	
POWER SUPPLY	100-240V; 50-60Hz Rechargeable batteries 3x NiMH AA1400mAh 1.2V	

Objective chart

Type	Magnification	Size (LxD)(mm)	Numerical aperture (N.A.)	Working distance (mm)	Coverglass thickness (mm)
Plan achromatic objectives	4x	28x19	0.1	37.5	0.17
	10x	43x20	0.25	6.54	0.17
	40x (R)	48.25x20	0.65	0.63	0.17
	100x (R)(I)	49x21.5	1.25	0.195	0.17

4. INSTALLATION

4.1 Instrument installation

- Take off the microscope from the carton box and remove all protections. Hold the microscope with both hands from the stand to avoid damaging any component and place it carefully on a stable flat surface. Remove all plastic and paper protections from each part of the microscope (eyepieces, light, diaphragm...)
- Place the head in its position and fix it with the lateral screw.
- Familiarize with the mechanical components of the microscope. Move these components to see how they behave and if they operate properly.