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Parameter list-2:

Parameter indicator	Name	Instruction for the parameter's function	(Setting range) Factory set value
Lc-	Password	When Lc=9, then we can see and modify parameters	0
Co-	Turn off the heat output deviation	When "PV ≥ SP + Co", turn off the heating output	(0.0 ~ 50.0 °C) 5.0
Hn-	Constant temperature time mode	0: minutes time; 1: hours time	(0 ~ 1) 0
En-	End of operation temperature	En = 0 end of run off output En = 1 end run to constant temperature	(0 ~ 1) 1
rH-	Range of temp. setting	The value of temperature setting	(0 ~ 100.0 °C) 70.0
SPL-	Lower limit	Temperature set value minimum value	(0 to highest limit) 0
SPH-	Highest limit	Temperature set value maximum value	((lower limit to highest limit) 70.0

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Environmental protection instructions

Do not dispose of this equipment with normal household waste at the end of its life cycle; take it to a collection point for the recycling of electrical and electronic equipment.

It does not contain elements that are hazardous or toxic to humans but improper disposal will harm the environment.

The materials are recyclable as indicated by the marking. By recycling materials or other forms of reuse of old equipment, you are making an important contribution to the protection of the environment. Please contact your local authority for advice on the collection points.

AT function

When the effect of the temperature control is not ideal, the AT function allows the system to be adjusted. Temperature auto tuning process can overshoot the set point; the users should be aware of this. In not running state, the controller will enter the auto-tuning of PID by pressing the “◀” button for 6s, “AT” indicator lights up, it will be not bright when the auto-tuning of PID is completed. When running the auto-tuning of PID, it can be stopped by pressing the “◀” button for 6s again. In the system auto tuning process, regardless of whether there is a constant temperature time setting, controller display window lower always displays the temperature setting value.

Internal parameters settings

Press the “SET” button for 3 seconds, controller will display the password prompt “Lc”. Adjust the password to the required value, then press the “SET” button again, it will run into the internal parameter setting state. If press the “SET” button for another 3 seconds, it will return to the running state.

Parameter list-1:

Parameter indicator	Name	Instruction for the parameter's function	(Setting range) Factory set value
Lc-	Password	When Lc=3 ,then we can see and modify parameters	0
AL-	Alarming setting	When the measured temperature is beyond “SP+AL”, the ALM indicator turns on. The buzzer sounds and the heater output turns off.	(0~100 °C) 5
T-	Control cycle	The heat control cycle of temperature	(1~60S) 5S
P-	Proportional band	Adjustment of proportional parameter	(1.0~rH) 26.5
I-	Integration time	Adjustment of integration parameter	(1~1000S) 415
d-	Differential time	Adjustment of differential parameter	(0~1000S) 415
Pb-	Zero point adjust	When the zero error is comparatively larger, this value should be updated. Pb=measured value –actual value	(-12.0~12.0 °C) 0.0
PK-	Full point adjust	When the full point error is comparatively larger, this value should be updated. PK=1000×(measured value –actual value)/ actual value	(-999~999) 0
Et-	Timing function	When Et = 0, no timing function; Et = 1 electric start timing; Et = 2 to the value set timing.	(0~2) 2

INDEX OF LANGUAGES

Spanish	2-7
English	8-12
French	13-18

INDEX OF CONTENTS

1. Safety warnings	8	4. Technical parameters	9
2. Product overview	9	5. Control panel.....	10
3. Product installation	9		

1. SAFETY WARNINGS

Prohibitions

- 1 It is not allowed to place the equipment in a place that is wet or may be splashed in water, which may cause short circuit or electric shock.
- 2 It is not allowed to insert the metal objects such as wires to any orifice and clearance of equipment, to prevent the occurrence of electric shocks or injuries.
- 3 It is not allowed to disassemble, repair, or modify the equipment by non-professionals.
- 4 Do not damage the power plug or power line. If the plug is loose, the power line must be replaced.
- 5 Do not place flammable or explosive materials near the incubator.

Observing items

- 1 The equipment should be placed on a firm surface or table. If the desktop is not strong enough or the location is not suitable, it will be caused injured.
- 2 The device must be reliably grounded in case of accidental leakage of electricity.
- 3 Be sure to turn off the power before any repairing or maintenance of this equipment, to avoid shock or injury.
- 4 Wearing gloves when making equipment maintenance, to avoid injury with the edges or sharp corners.
- 5 Do not damage the power cord or use a non-specified power cord.
- 6 If the device is found to be abnormal, immediately remove the plug and terminate the device. It may cause shock or fire in abnormal conditions.

Attention items

1. The equipment should be placed on a level surface.
2. Use an independent power socket with ground wire and insert the plug tightly.
3. Unplug the power cord before moving the incubator to a new location.
4. Internal parameters setting should be responsible by the specific management personnel, in order to prevent the controller program dysfunction.
5. Equipment installation location should have good ventilation. Should be placed away from the wall; the distance from objects must be greater than 30 cm.
6. The surfaces shall be not contacted with the volatile chemicals such as gasoline, thinner, etc.
7. Keep the incubator inside and outside clean.

2. PRODUCT OVERVIEW

This device is suitable for microorganism culture; it is an important equipment used in medicine, biology, biochemistry, etc.

Operating principle

This incubator uses an electric resistance for heating. It has a special design that ensures a very good temperature uniformity inside the working chamber and an accurate temperature control.

Main characteristics

1. The shell of this incubator is made of high quality cold-rolled steel plate. The surface is protected with a corrosion resistant and durable coating.
2. The working chamber is made of mirror stainless steel, easy to clean.
3. Double door structure: high quality tempered glass and outer door, both with rubber seals.
4. Microcomputer based intelligent control system, PID temperature control, LED screen, 4 buttons control panel, with over-temperature alarm and timing function.


3. PRODUCT INSTALLATION

1. The incubator should be placed far away from the dust and vibration place.
2. Ensure the four rubber legs are evenly on the ground.
3. The shell must have a reliable grounding (there is a line for the grounding in the three-pin plug).
4. After installing the incubator, open the doors to take out the perforated trays and insert them in the appropriate level according to the experiment requirements.

4. TECHNICAL PARAMETER

Code	JBB006	JBB007	JBB008
Temperature range	RT+5~65 °C		
Temperature resolution	0,1 °C		
Temperature fluctuation	±0,5 °C		
Temperature uniformity	±1 °C		
Chamber material	Mirror stainless steel		
Shell material	Coated cold-rolled steel plate		
Insulation material	Polyurethane (PU)		
Heater	Heating resistance, tube form		
Rated power	0.25kw	0.5kw	0.6kw
Air vent	Inner diameter 28mm*1, located at the top		
Temperature displays			
Temperature display	4 digits upper display (PV)		
Set temperature display	4 digits lower display (SV)		
Timer	0001~9999 minutes (or hours) / Continuous mode 0000		
Temperature sensor	Pt100		
Safety device	Over-temperature alarm		
Inner size (LxAxH mm)	350x350x350	400x350x450	500x450x550
Exterior size (LxAxH mm)	500x475x615	550x475x715	650x575x815
Capacity	43L	63L	124L
Max. shelf load	15 kg		
Shelves included	2		
Voltage (50/60Hz) / Rated current	AC220V/1.1A	AC220V/2.3A	AC220V/2.7A
Net weight, kg	29	33	49

5. CONTROL PANEL

Light indicators

- 1) "AT" indicator: When the controller enters the auto-tuning of PID, this indicator is bright.
- 2) "OUT" indicator: If the heater output turns on, this indicator is bright, else this indicator is not bright.
- 3) "ALM" indicator: When the over-temperature alarm occurs, this indicator is bright.


Operation

1) When the controller is switched on, display windows show "InP" and the value of temperature range for 3 seconds, then it starts running.

2) Temperature and time settings:

Press the "SET" button, the controller runs into the temperature setting state. Re-press the "SET" button, the controller runs into the time setting state. In setting state, you can use the "◀", "▼" and "▲" buttons to get the required settings. Press the "SET" button again, it returns from the setting state and the settings are saved automatically.

If the time is set as "0", the incubator will run continuously, the display window of "SV" will display the set point temperature. If the time set value is not equal "0", timer starts when the measured temperature reaches the set point temperature, the display window of "SV" will display the runtime. When the runtime is over, the "SV" window will display "End", the buzzer will sound for 30s; press the button "▼" for 3s, the program will restart.

(1) The normal display

(2) The temperature setting state

(3) The time setting state

(4) Timing display state

3) When temperature alarm, the buzzer will sound, "ALM" lights. If a change in temperature setting and over-temperature alarm, "ALM" lights up, but no songs buzzer.

4) When the buzzer sounds, it can be muted by pressing any button.

5) "◀" button: In the setting state, it changes the digit to be adjusted.

6) "▼" button: In the setting state, it can reduce the set value by pressing the button. If press and hold the button, the set value will reduce continuously. The timing state, long press the button for 3 seconds and can make the program stop.

7) "▲" button: In the setting status, it can increase the set value by pressing the button. If press and hold the button, the set value will increase continuously.

8) In setting state, the controller will return to run status if without any key press in one minute.

9) If the display window shows "----", it indicates the fault of temperature.