**AGITADOR MAGNÉTICO DIGITAL 20 L
DIGITAL MAGNETIC STIRRER 20 L
AGITATEUR MAGNETIQUE NUMERIQUE 20 L**

Modelo / Model / Modèle / RSLAB-4C

RS Lab

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INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, does not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non-adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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INDEX OF CONTENTS

1. Safety instructions	14	8. Function: stirring	20
2. Proper use	15	9. Remote control	20
3. Inspection	15	10. Faults	20
4. Trial run	16	11. Maintenance and Cleaning	21
5. Control and Display	16	12. Associated Standards and regulation ..	21
a. Control elements	16	13. Technical data	22
b. Display	17	14. Accessories	23
6. Initialization	18		
7. Function: heating	19		
a. External temperature sensor	19		
b. Residual heat warning (Hot)	19		
c. Safety temperature	19		

Welcome to the digital magnetic stirrer RSLAB-4C (heating) user manual. Users should read this manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument.

When help needed, you can always contact the service department of manufacturer or your supplier for technical support. Please provide the customer care representative with the following information: Serial number, description of problem, methods and procedures adopted to resolve the problems and your contact information.

These RSLAB stirrers are warranted to be free from defects in materials and workmanship under normal use and service for a period of 24 months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation.

For claims under the warranty please contact your local dealer. You may also send the instrument direct to our works, enclosing the invoice copy and giving reasons for the claim.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

1- SAFETY INSTRUCTIONS

- ◆ Read the operating instructions carefully before use the instrument.
- ◆ Ensure that only trained staff work with the instrument.

Risk of burn

- ◆ Caution when touching the housing parts and the heating plate. The heating plate can reach temperatures of 550 °C.
- ◆ Pay attention to the residual heat after switching off.

Protective ground contact

- ◆ Make sure that socket is earthed (protective ground contact) before use.
- ◆ When work, wear the personal guard to avoid the risk from splashing and evaporation of liquids and release of toxic or combustible gases.
- ◆ Set up the instrument in a spacious area on a stable, clean, non-slip, dry and fireproof surface, do not operate the instrument in explosive atmospheres, with hazardous substances or under water.
- ◆ Gradually increase the speed, reduce the speed if the stirring bar breakaway because of too high speed and if the instrument is not running smoothly, or if the container moves on the base plate.
- ◆ Temperature must always be set to at least 50 °C lower than the fire point of the media used.
- ◆ Beware of hazards due to flammable materials or media with a low boiling temperature, overfilling of media and unsafe container.
- ◆ Process pathogenic materials only in closed vessels.
- ◆ Check the instrument and accessories before handle it for damaged component. Safe operation is only guaranteed with the accessories described in the "Accessories" chapter. Accessories must be securely attached to the device and cannot come off by themselves. Always disconnect the plug before fitting accessories.
- ◆ Ensure that the external temperature sensor is inserted in the media to a depth of at least 20 mm.
- ◆ When using metal vessels, do not place the temperature sensors on the bottom of the vessel. Placing sensors on the vessel bottom can cause excessively high

14- ACCESSORIES

[1] 90694001

[2] 90694002

[3] 90694011

[4] 90694305

[5] 90694306

[6] 90694307

[7] 90694308

[8] 90694309

[9] 90694310

Reference

- [1] 90694001
- [2] 90694002
- [3] 90694011
- [4] 90694305
- [5] 90694306
- [6] 90694307
- [7] 90694308
- [8] 90694309
- [9] 90694310

Description


- Temperature sensor pt1000, l. 230 mm, **included**
- Stand kit for sensor pt1000, **included**
- Temperature sensor pt1000, with glass coated, l. 230 mm
- Poly-bloc stand red, 11 holes \varnothing 12,5 mm / width 20 mm / tubes 4 ml
- Poly-bloc stand purple, 4 holes \varnothing 28 mm / width 24 mm / tubos de 20 ml
- Poly-bloc stand blue, 4 holes \varnothing 28 mm / width 30 mm / tubos de 30 ml
- Poly-bloc stand black, 4 holes \varnothing 28 mm / width 43 mm / tubos de 40 ml
- Poly-bloc stand green, 6 holes \varnothing 17,8 mm / width 26 mm / tubos de 8 ml
- Poly-bloc stand gold, 4 holes \varnothing 21,6 mm / prof. 31,7 mm / tubos de 16 ml

13- TECHNICAL DATA

<i>Voltage</i>	200-240/100-120 VAC
<i>Frequency</i>	50/60 Hz
<i>Power</i>	1050 W (heating) / 50 W
<i>Stirring point position quantity</i>	1
<i>Max. stirring quantity (H2O)</i>	20 Liters
<i>Max. magnetic bar</i>	80x10 mm
<i>Motor</i>	Brushless motor
<i>Max. power input of motor</i>	18 W
<i>Max. power output of motor</i>	10 W
<i>Speed range</i>	100-1500 rpm
<i>Rotary speed display</i>	LCD
<i>Speed display accuracy</i>	1 rpm
<i>Hotplate material</i>	Glass ceramic
<i>Dimensions of workplate</i>	184x184 mm
<i>Dimensions of heating</i>	176x176 mm
<i>Heating power</i>	1000 W
<i>Temperature range</i>	RT-550 °C
<i>Temperature display</i>	LCD
<i>Temperature display accuracy</i>	0,1 °C
<i>Control accuracy of heating temperature</i>	1 °C
<i>The safety temperature range of the hotplate</i>	100-580 °C
<i>Temperature sensor in medium</i>	PT1000
<i>Control accuracy of heating temperature sensor</i>	0,2 °C
<i>Residual heat warning</i>	50 °C
<i>Dimensions</i>	215x360x112 mm
<i>Weight</i>	5,3 kg
<i>Permitted ambient temperature</i>	5-40 °C
<i>Permitted relative humidity</i>	80 %
<i>Protection class acc. to DIN 60529</i>	IP21
<i>RS232 Interface</i>	Yes

temperatures to be measured especially in media which have poor conductivity. The tip of the measuring sensor must be at least 5 mm from the vessel bottom; a distance of 10 mm is ideal.

- ◆ The instrument can only be disconnected from the main power supply by pulling out the main plug or the connector plug.
- ◆ The voltage started on the label must correspond to the main power supply.
- ◆ Ensure that the main power supply cable does not touch the heating base plate. Do not cover the device.
- ◆ The instrument may only be opened by experts.
- ◆ Keep away from high magnetic field.
- ◆ Observe the minimum distances between the devices, between the instruments and the wall and above the assembly (min. 100 mm).

**2- PROPER USE**

The instrument is designed for mixing and/or heating liquids in schools, laboratories or factories. This device is not suitable for using residential areas or other constraints mentioned in chapter 1.

3- INSPECTION

- ◆ Unpack the equipment carefully and check for any damages which may have arisen during transport. If it happens, please contact manufacturer for technical support.
- ◆ If there is any apparent damage to the system, please do not plug it into the power line.

◆ The RSLAB-4NC y RSLAB-4C includes the following items:

- Main unit	1 unit	- Stirrer bar	1 unit
- Temp. sensor PT1000	1 unit	- Power cable	1 unit
- Kit for sensor PT1000	1 unit	- User manual	1 unit


4- TRIAL RUN

- ◆ Make sure the required operating voltage and power supply voltage match.
- ◆ Ensure the socket must be earthed reliably.
- ◆ Ensure the power is off.
- ◆ Plug in the power cable, ensure the power be on and begin initializing.
- ◆ Add the medium into the vessel with a stirring bar.
- ◆ Put the vessel on the work plate.
- ◆ Set the rated stirring speed and begin.
- ◆ Observe the stirring bar and LCD display.
- ◆ Set the rated temperature and start heating.
- ◆ Observe the real temperature on LCD display.
- ◆ Stop the heating and stirring functions.
- ◆ If these operations above are normal, the device is ready to operate. If these operations are not normal, the device may be damaged during transportation, please contact manufacturer for technical support.

5- CONTROL AND DISPLAY

a) Control elements

1. LCD display
2. LED heat display
3. LED stir display
4. Heating knob
5. Stirring knob
6. Switch ON/OFF
7. Safe Temperature Button

11- MAINTENANCE AND CLEANING

Warning! Cut off power when maintenance and cleaning

- ◆ Proper maintenance can keep instruments working in a good state and lengthen its lifetime.
- ◆ Do not spray cleanser into the instrument when cleaning.
- ◆ Unplug the power line when cleaning.
- ◆ Only use recommended cleaners:
 - Dyes: isopropyl alcohol
 - Construction materials: water containing tenside / Isopropyl alcohol
 - Cosmetics: water containing tenside / Isopropyl alcohol
 - Foodstuffs: water containing tenside
 - Fuels: water containing tenside
- ◆ Wear proper protective gloves during cleaning of the instrument.
- ◆ Before using other method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does no destroy the instrument.
- ◆ Unit should be clean and fitted in the original packing or a suitable packing for transport.
- ◆ The enamel makes the heating plate easier to care for and more resistant to acids and bases. Because of it, however, the heating plate is also more susceptible to extreme fluctuations in temperature and the force of impact. This can result in cracks forming or the coating flaking off.

12- ASSOCIATE STANDARDS AND REGULATIONS

Construction in accordance with the following safety standards:

- EN 61010-1
- UL 3101-1
- CAN/CSA C22.2 (1010-1)
- EN 61010-2-10

Construction in accordance with the following EMC standards:

- EN 61326-1

Associated EU guidelines:

- EMC guidelines: 89/336/EWG
- Instrument guidelines: 73/023/EWG

8- FUNCTION: STIRRING

- ◆ The function “**stirring**” is switched on or off via the rotating stirring knob. The motor speed is set on the knob (100 to 1500 rpm in steps of 10 rpm).
- ◆ When both functions heating and stirring are switched on and those above operations are done, the LCD will shift to the speed value and come back to the temperature value in about 5 seconds.

9- REMOTE CONTROL

The unit can be controlled from an external PC (using the dedicated software) via the RS 232 C serial interface fitted to the unit. Data communication from laboratory instrument to computer is only possible on demand of the computer.

10- FAULTS

- ◆ Device can't be power ON:
 - Check whether the power cable is plugged
 - Check whether the fuse is broken or loose
- ◆ Fault in power on self-test:
 - Switch OFF the unit, then switch on and reset the instruments to factory default setting
- ◆ Temperature cannot reach set point:
 - Check whether the safety temperature value is set too low
- ◆ Stir speed cannot reach the set point:
 - Check whether the stirring bar breakaway monitoring is ON
 - Excessive medium viscosity may cause abnormal speed reduction of the motor
- ◆ Heating or stirring cannot be started when push the temperature control knob:
 - Check whether the unit is in the mode setting program
 - Exit from the mode setting and restart the heating/stirring functions
- ◆ Unit cannot be powered off when switched off:
 - Check if the residual heat warning functions is still ON and hotplate temperature is above 50 °C (The LCD still work and “**Hot**” flash)
 - Turn off residual heat warning function before powered OFF the unit

If these faults are not resolved, please contact manufacturer for technical support.

Items	Description
Speed control Knob	Set the rated rotary speed from 0 to 1500 rpm. The stirring function is switched on or off by pushing the knob
Temp. control Knob	Set the rated temperature. The heating function is switched on or off by pushing the knob
LCD Display	The LCD screen shows the actual working state and all selected parameters
LED Heat	When the heating function is switched on, the LED heat is lit
LED Stir	When the stirring function is switched on, the LED stir is lit
Safe Temp. Button	Set safety temperature, stop heating when exceed secure temperature
Power Switch	Switch on or off the instrument

ENGLISH

b) Display

Character	Description
Temp and C°	Display temperature when the heating function is switched on
Mot	Display stirring state when the stirring function is switched on
Hot	Display hot warning if the hotplate temperature is above 50 °C after switching off the heating function
Probe	Display when using external probe
Remote	Display in case of remote control
Err	Display in case of error happening
Rated value Real value	Display value in case of heating or stirring function switching on

If both heating and stirring functions have been started at the same time, heating function always has higher priority. If in this case speed is changed via the stirring knob, it displays stirring speed and reverses to temperature in the duration of 5 seconds.

6- INITIALIZATION

- ◆ Turn on instrument by the power switch.
- ◆ The “heat” and “stir” LEDs and the LCD display flash about 3 times.
- ◆ LCD displays the workplate safe temperature limit in “SAFE xxx °C”.
- ◆ LCD displays residual heat warning “rES On/OFF”
- ◆ LCD displays stirring bar breakaway monitoring “br On/OFF”
- ◆ LCD displays the “Probe” characters in case of connecting with external temperature sensors.
- ◆ RSLAB-4C are supplied with the following default settings:
 - Setting temperature: 25 °C
 - Setting rotary speed: 100 rpm

7- FUNCTION: HEATING

- ◆ The device is controlled by digital temperature control technology, which has two separate safe circuits. The hotplate is kept at a constant temperature by a digital control circuit. The hotplate temperature can also be monitored from a separate, adjustable safe circuit. The two temperature sensors (Pt1000) internal for temperature control are built into the hotplate. The single external Pt1000 can monitor the temperature of sample.
- ◆ If external temperature sensors needed, it must be plugged in before the device is switched on.
- ◆ Ensure the external temperature sensor connect firmly before heating.
- ◆ Set the temperature via turning the temperature control knob slowly to the ideal value.
- ◆ When the heating function is switched in, the LED “heat” will light up and the LCD will display the real temperature.
- ◆ The rated temperature will be displayed on the right-hand side of the LCD as well as Temp and °C characters.
- ◆ The heating function is switched on or off by pushing heating knob.

a. Working with external temperature sensor

- ◆ The external temperature sensor PT1000 is a standard accessory. If the external temperature sensor is not plugged in, “Probe” will flash. If plugged in, “Probe” will be shown on the digital display to indicate that the sensor is operating. The external temperature sensors’ value and actual temperature are displayed. Compared with the temperature control of the hotplate, the external temperature sensor can control the media temperature precisely.
- ◆ Do not insert or remove sensor while heating.

b. Residual heat warning (HOT)

- ◆ In order to prevent the risk of burns from a hotplate, RSLAB-4C has a residual heat warning function. When the heating function is switched off and the hotplate temperature is still above 50 °C, “hot” will flash to warn that there is a hazard of burns from the hotplate.
- ◆ When the heating plate temperature drops to below 50 °C, the unit will automatically switch off. When the plug is pulled out, the residual heat warning function cannot be run.

c. Setting the safety temperature

- ◆ The safety temperature can be set adjusting the temperature setting key to set the safety temperature.