**VASOS COMUNICANTES
COMMUNICATING VESSELS**

REF : QLB004

INDEX OF LANGUAGES-ÍNDICE DE IDIOMAS

ESPAÑOL2
ENGLISH.....3

VASOS COMUNICANTES

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COMUNICATING VESSELS

This instrument is composed by four glass tubes with different diameter and shape. The tubes are connected on the base through a horizontal glass tube that has a vertical one in the middle that is sealed with a platform. This instrument allow to demonstrate that the level of liquid is the same in all tubes it does not depend on the diameter and shape.

Recommendations of use for school

Is recommended for secondary school especially for 4º and 2º. The theory that they can learn is the principle of Pascal "The pressure exerted in a liquid point is transmitted with the same intensity in all directions".

The applications can be communicant vessels and hydraulic press.

In this case we describe two tips for the use of communicant vessels at school. First, the students can prove that the level of the liquid is the same in all vessels, they can choose different liquids and describe what happened.

In addition to see the instrument at the laboratory they can search about the real use of this principle for example the location of water reservoir respect the town, is the same structure, in this case the pipes represent the communicant vessels due to the water reservoir must be in a high level.