PLANO INCLINADO

INCLINED PLANE

REF. QLB001

ÍNDICE DE CONTENIDO – INDEX OF LANGUAGES

ESPAÑOL.....2

ENGLISH.....4

PLANO INCLINADO

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INCLINED PLANE

The instrument has been designed to study concepts of motion on an inclined Surface and its relationship to the friction. All metal construction provides sturdiness and durability to this apparatus. It comprises of two metal boards, hinged at one end. One board is used as base another one as inclined plane. One side of inclined board has a linear scale graduated at every 5 cm to measure linear distance along its length and it can be inclined and fixed at any angle from 0-40 to the horizontal measured on the attached angle measurer. A vertical scale 0-25 cm is also printed on the angular scale to measure the height of inclined surface at any point. The end of the inclined board features an adjustable pulley wheel from which masses, weight pan etc. can be suspended. A plastic trolley with hook and a pan with lightweight cord are also provided with the incline plane.

Object: to measure the coefficient of friction, making use of an inclined plane; and to study equilibrium and non-equilibrium of a body on an inclined plane under the reaction of forces.

Apparatus required: inclined plane, plastic trolley, weight pan hanger, set of weights, clock timer and meter stick.

Theory:

Sliding friction: it is the friction between two surfaces of the bodies in sliding motion.

Force of sliding friction: it is the least force required to make a body start sliding over a surface.

$$F \propto R$$

$$F = \mu R$$

Force friction

Where μ = coefficient of friction, R is the normal reaction.

At equilibrium

$$F = P + p$$

$$R = W + w$$

Procedure

1. Weigh the plastic trolley and put it on the table top
2. Tie one end of a thread with the hook of the plastic trolley and pass it over the pulley.
3. Find the weight of the pan.
4. Tie other free end of the thread with the pan and let the pan hang vertical. The pan will pull the plastic trolley horizontally by a force equal to its weight.
5. Since the pan itself does not pull the plastic trolley, put some weight in the pan.

6. Increase weights in the pan little by little, till the plastic trolley just starts rolling on the table top.
7. Note the total weights put in the pan then record those in observation table (sum of weight of pan and weights in pan gives the force of sliding friction).
8. Put one 50g or 20g weight over the pan and repeat steps 6 and 7.
9. Repeat steps 6, 7, 8 six times. Every time increase weight by 50g or 20g.
10. Also repeat the experiment by placing same weight in the plastic trolley
11. Record the observations in tabular form as given below.

Observations:

Weight of plastic trolley	$W = \dots \text{g-wt}$
Weight of pan	$P = \dots \text{g-wt}$

Calculations: total weight (force) pulling the plastic trolley and weights gives the value of force of sliding friction.

On horizontal surface, total weights being pulled give normal reaction R.

Total weight (force) pulling these weights gives dynamic friction F.

Draw a graph between normal reaction R and limiting friction F taking R along X axis and F along Y axis. The graph comes out to be straight line.

Results: It is found that as the total weight pulled increases, force of the limiting friction also increases. This increase is in the direct proportion.

The graph shows that the limiting friction F is directly proportional to the normal reaction R. it is an agreement with the law of limiting friction.

The constant ratio F/R is called coefficient of friction. It can be calculated by finding the slope of the graph.

Precautions:

1. The surface should be horizontal
2. The thread part between plastic trolley and pulley should be horizontal
3. Weight in the pan should be increased in small steps and pan should not oscillate or rotate.
4. Table top should be tapped gently each time.
5. Pulley should be friction less.