



## Squibb with teflon key

Code PQB001, PQB002, PQB003, PQB004, PQB005, PQB006, PQB007, PQB008, PQB009, PQB010

- Separatory funnels are used for the separation of immiscible liquids of different densities, typically water and various organic solvents (ether, chloroform, or alcohols).
- Made of borosilicate glass resistant to thermal shock, Endo Glassware's separatory funnels feature a pear shape and a plastic stopper.
- The stopcock made of Teflon is resistant to the action of acids and bases.



Code	PQB001	PQB002	PQB003	PQB004	PQB005	PQB006	PQB007	PQB008	PQB009	PQB010
Capacity (mL)	50	100	125	250	250	250	500	500	1000	2000
Ø Funnel (mm)	55	65	65	80	80	80	100	100	125	150
Ø Branch (mm)	8	8	8	8	8	8	10	10	10	12
Height (mm)	230	245	250	300	300	300	350	350	400	525
Grinding	19/26	19/26	19/26	29/32	19/26	24/29	29/32	24/29	29/32	29/32